

**Shosaikoto Extract**

小柴胡湯エキス

**Change the Identification (1) as follows:**

**Identification (1)** Shake 2.0 g of dry extract (or 6.0 g of the viscous extract) with 10 mL of sodium hydroxide TS, then add 5 mL of 1-butanol, shake, centrifuge, and use the supernatant liquid as the sample solution. Separately, dissolve 1 mg of saikosaponin b<sub>2</sub> for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot 10  $\mu$ L of the sample solution and 2  $\mu$ L of the standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, ethanol (99.5) and water (8:2:1) to a distance of about 10 cm, and air-dry the plate. Spray evenly 4-dimethylaminobenzaldehyde TS for spraying on the plate, heat at 105°C for 5 minutes, and examine under ultraviolet light (main wavelength: 365 nm): one of the spot among the several spots obtained from the sample solution has the same color tone and R<sub>f</sub> value with the yellow fluorescent spot obtained from the standard solution (Bupleurum Root).

*Suppl II*  
*JPXVI*  
*(2014)*

*Suppl II*  
*JPXVI*  
*(2014)*

*Suppl II*  
*JPXVI*  
*(2014)*

*Suppl II*  
*JPXVI*  
*(2014)*

*Suppl II*  
*JPXVI*  
*(2014)*