Suppl I, JP XVIII (2022)

Suppl I, JP XVIII (2022)

Suppl I, JP XVIII (2022)

## Keishibukuryogan Extract

桂枝茯苓丸エキス

## Change the Assay (3) as follows:

Assav

(3) Amygdalin—Weigh accurately about  $0.5 \, \mathrm{g}$  of the dry extract (or an amount of the viscous extract, equivalent to about  $0.5 \, \mathrm{g}$  of the dried substance), add exactly  $50 \, \mathrm{mL}$  of diluted methanol (1 in 2), shake for 15 minutes, filter, and use the filtrate as the sample solution. Separately, weigh accurately about  $10 \, \mathrm{mg}$  of amygdalin for assay, dissolve in diluted methanol (1 in 2) to make exactly  $50 \, \mathrm{mL}$ , and use this solution as the standard solution. Perform the test with exactly  $10 \, \mu \mathrm{L}$  each of the sample solution and standard solution as directed under Liquid Chromatography  $\langle 2.01 \rangle$  according to the following conditions, and determine the peak areas,  $A_{\mathrm{T}}$  and  $A_{\mathrm{S}}$ , of amygdalin in each solution.

Amount (mg) of amygdalin =  $M_S \times A_T/A_S$ 

 $M_{\rm S}$ : Amount (mg) of amygdalin for assay taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 210 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silicagel for liquid chromatography (5  $\mu$ m in particle diameter).

Column temperature: A constant temperature of about  $45\,^{\circ}\mathrm{C}$ .

Mobile phase: A mixture of 0.05 mol/L sodium dihydrogen phosphate TS and methanol (5:1).

Flow rate: 0.8 mL per minute (the retention time of amygdalin is about 12 minutes).

System suitability—

System performance: When the procedure is run with  $10 \mu L$  of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of amygdalin are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with  $10 \,\mu\text{L}$  of the standard solution under the above operating conditions, the relative standard deviation of the peak area of amygdalin is not more than 1.5%.